



Instruction for Use

***Evo M-MLV* RT Kit for qPCR**

AG11707

Version.V4E1

**Research Use Only
Not For Diagnosis Procedures**

1. Description

This product is a reverse transcription reagent specifically designed for Real-Time RT-PCR. It incorporates the high-processivity Evo M-MLV reverse transcriptase, enabling rapid and efficient cDNA synthesis in a short time. The resulting cDNA is fully compatible with both dye-based and probe-based qPCR assays, and can be used with appropriate downstream reagents to achieve high-performance gene expression analysis tailored to different experimental objectives.

2. Kit Information

Kit Name	Cat. No	Specification
<i>Evo M-MLV</i> RT Kit for qPCR	AG 11707	200 rxns / 10 μ l

3. Transportation and Storage

Storage	Store at -20°C
Transportation	Transport at -20°C Dry Ice or Blue Ice Condition

4. Kit Components

Kit Components	Volume
<i>Evo M-MLV</i> RTase Enzyme Mix* ¹	100 μ l
5X RTase Reaction Buffer Mix II* ²	400 μ l
Oligo dT (18T) Primer (50 μ M)	100 μ l
Random 6 mers Primer (100 μ M)	400 μ l
RNase Free Water	1 ml x 2 Pcs

*1: Contains *Evo M-MLV* RTase, RNase Inhibitor.

*2: Contains dNTP.

5. Precautions

- 1) Before use, briefly centrifuge the *Evo M-MLV* RTase Enzyme Mix to collect all enzyme solution at the bottom of the tube and prevent loss of enzyme during handling.
- 2) When setting up multiple reactions, first prepare a master mix of all required reagents, then aliquot it into individual reaction tubes.
- 3) All reaction mixtures should be prepared on ice.

6. Protocol

Prepare the RT mix as below

Components	Input Volume*2
5X RTase Reaction Buffer Mix II	2 μ l
Evo M-MLV RTase Enzyme Mix	0.5 μ l
Oligo dT (18T) Primer (50 μ M)*1	0.5 μ l
Random 6 mers Primer (100 μ M)*1	0.5 μ l
Total RNA*3	-
RNase free water	Up to 10 μ l

Reaction Program:

Temperature	Time
37°C*4	15 min
85°C	5 sec
4°C	Hold*5

- When using the SYBR Green qPCR method, add 0.5 μ l Random 6 mers Primer (100 μ M).
 When using the probe-based qPCR method, add 2 μ l Random 6 mers Primer (100 μ M).
 Alternatively, the primer mix can be omitted and individual primers may be used as follows:
 Oligo dT (18T) Primer: 25 pmol per 10 μ l reaction
 Random 6 mers Primer: 50 pmol per 10 μ l reaction (for SYBR Green qPCR)
 Random 6 mers Primer: 200 pmol per 10 μ l reaction (for probe-based qPCR)
 Gene-Specific Primer: 2.5 pmol per 10 μ l reaction (adjustable within 1-5 pmol)
- To ensure accurate preparation of the reverse transcription reaction, it is recommended to prepare a master mix of all components first, aliquot into reaction tubes, and finally add the RNA samples. The reaction volume may be adjusted as required.
- The amount of RNA can be adjusted as needed. In a 10 μ l reverse transcription reaction, use up to 500 ng total RNA for SYBR Green qPCR and up to 1 μ g total RNA for probe-based qPCR.
- When using a Gene-Specific Primer, the reaction may be performed at 42 °C for 15 min. To reduce non-specific amplification, the reaction temperature can be increased to 50°C.
- 5: If the reaction product is used immediately for subsequent qPCR, it may be kept at 4°C or on ice. For short-term storage, store at -20°C; for long-term storage, store at -80°C.

7. Quantitative PCR Reaction Analysis

The reaction mixture obtained above can be directly used for subsequent quantitative PCR. Ensure that the added volume does not exceed 1/10 (V/V) of the qPCR reaction volume.



Accurate Biotechnology (Hunan) Co., Ltd

No.1 JinYang Road,
WangCheng District, Changsha City, Hunan Province, China
service@agbio.com.cn
+86 400 767 6022
en.agbio.com.cn

Research Use Only